




Deflector for adapter 61808250

6.2058.070

The deflector is mounted with the fastening plate 6.2058.050 at the side of the Sample Processor tower and allows the automatic pulling off the used pipetting tips into the collection container 6.1625.010.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery 6.2058.070

Qt.	Order no.	Description
1 PCS	6.2621.100	3 mm hex key for IC Sample Processors



Optional accessories

Order no.	Description
6.1808.250	Pipette adapter 1-10 mL for transfer head Adapter with outside thread M6 and Eppendorf conus to use with transfer head 6.1462.240 for usage of 10 mL Eppendorf tips

