



## Liquid Handling Station, right

6.2841.130

Combined rinsing and dilution station for sample changers with two towers and a Swing Head. To be assembled to the right of the tower.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

### Scope of delivery 6.2841.130

| Qt.   | Order no.  | Description        |
|-------|------------|--------------------|
| 1 PCS | 6.2621.140 | Hexagon key 2.5 mm |



1 PCS

6.1808.050

Tubing adapter olive / M8 outer

1 M8 outer thread and 1 tubing olive. E.g. for thermostat jacket of exchange units and stability measuring instruments.



---

1 PCS

6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m



---

1 PCS

6.2621.030

Hexagon key 4 mm



---

2 PCS


6.2744.060

Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



# Optional accessories

| Order no.  | Description                                                                                                                                                     |                                                                                       |
|------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| 6.2762.100 | Mixing vessel for Liquid Handling Station<br><br>Replacement mixing vessel for combined rinsing and dilution stations.                                          |    |
| 6.2762.200 | Transparent cover for the 6.2762.100 mixing vessel, PMMA<br><br>Transparent cover for the 6.2762.100 mixing vessel, for visual inspection in the mixing vessel. |    |
| 2.741.0010 | 741 Magnetic Stirrer<br><br>Magnetic stirrer for sample changers.<br><br>The magnetic stirrer is located beneath the sample rack.                               |  |