



Equipment for pH Adjustment at the 845 Eluent Synthesizer

6.5329.000

pH adjustment to Eluent Synthesizer. Includes pH electrode and additional accessories.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.5329.000

Qt.	Order no.	Description
1 PCS	6.1446.000	SGJ stopper / B-14/(15)






Eluent bottle head to connect up to 4 tubings with M6 connectors. For eluent bottle of the Eluent synthesizer.



Connection eluent bottle at the eluent Synthesizer - Titrino



Additional remote interface for the connection of devices that can be controlled via remote lines. With permanently attached cable.



1 PCS

6.2726.090

Antidiffusion valve

For antidiffusion buret tips 6.1541.XXX, 6.1543.130, 6.1543.140,
6.1543.150, 6.1543.160, 6.1543.200



1 PCS

6.3026.210

Exchange Unit 10 mL

Exchange unit with integrated data chip with 10 mL glass cylinder and light protection. PCTFE/PTFE flat stopcock, FEP tubing connection, antidiffusion buret tip and standard amber glass reagent bottle.

