



IC equipment : PCR with dosing unit

6.5330.400

The IC equipment: PCR with dosing unit makes it possible to convey a PCR reagent with the 800 Dosino, for example for the photometric determination of UV/VIS-active components after post-column derivatization.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.5330.400

Qt.	Order no.	Description
2 PCS	6.1608.030	Round glass bottle / 1000 mL / GL 45



3 PCS

6.1805.110

FEP tubing / M6 / 80 cm

With light and kink protection



2 PCS

6.1619.000

Adsorber tube for 807 Dosing Unit



1 PCS

6.2744.010

Pressure screw 5x

With UNF 10/32 connection. For the connection of PEEK capillaries



1 PCS

6.2744.080

M6 thread / UNF 10/32 adapter

With M6 outer thread. For IC instruments.



2 PCS

6.1819.100

FEP aspiration tubing to 2 L bottle

Aspiration tubing to bottle attachment with M6 connector. For the production of eluents.



2 PCS

6.1602.160

Bottle attachment GL 45 for eluents and reagents

For eluent, reagent, and rinsing bottles, with connectors for the adsorber tube and the aspiration tubing.



1 PCS

6.2744.200

Adapter UNF 10/32 outer / M6 inner

Connecting 1/16 in. PEEK capillaries and transfer tubing (M6)



2 PCS

6.1624.000

Adapter SGJ 14 for 6.1619.XXX Adsorber tube

Adapter SGJ 14. Part of 6.1619.010 for Dilution Sample Processor.



1 PCS

6.1831.030

PEEK capillary 0.75 mm i.d. / 3 m

For all IC components.



1 PCS

6.2057.250

Dosino holder for PCR

Dosino holder for IC instruments for PCR Dosino.



1 PCS

6.3033.150

Dosing unit PCR 5 mL

5 mL dosing unit specially suitable for conveying solutions for post-column derivatization (PCR) in IC.

