



# VA electrode equipment with Multi-Mode Electrode pro for Professional VA instruments

6.5339.030

Complete electrode set for polarographic and voltammetric determinations. Includes Multi-Mode Electrode pro, reference electrode, platinum auxiliary electrode, measuring vessel, stirrer, electrolyte solution and additional accessories for setting up and operating the Multi-Mode Electrode.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

## Scope of delivery 6.5339.030

| Qt. | Order no. | Description |
|-----|-----------|-------------|
|-----|-----------|-------------|

1 PCS

6.0343.100

Separate Pt Rod Electrode

Platinum auxiliary electrode for voltammetry and CVS.



---

1 PCS

6.0728.120

Ag/AgCl reference electrode

Inner reference system with electrolyte  $c(\text{KCl}) = 3 \text{ mol/L}$ . For use with electrolyte vessel made of plastic 6.1245.010.



---

1 PCS

6.1204.500

Stirrer for Professional VA instruments

Stirrer for Professional VA instruments, complete, metal-free.



---

1 PCS

6.1226.030

Mercury drop capillary (non-silanized) for MME

Non-silanized mercury drop capillary (10 pieces)



1 PCS

6.1226.050

Mercury drop capillary (silanized) for MME

Silanized mercury drop capillary (10 pieces)



1 PCS

6.1245.010

Electrolyte vessel made of plastic

Electrolyte vessel made of plastic with ceramic diaphragm for double junction reference electrode.




1 PCS

6.1246.120

Multi-Mode Electrode pro

Mercury electrode for voltammetry. Can be operated as DME, SMDE or HMDE.



1 PCS

6.1247.020

Needles for MME

Needles for Multi Mode Electrode (3 pieces).



1 PCS

6.1415.210

Titration vessel / 10-90 mL



---

1 PCS

6.2308.020

Electrolyte 3 mol/L KCl (250 mL)

Electrolyte solution  $c(\text{KCl}) = 3 \text{ mol/L}$ , (for Ag/AgCl reference systems)



---

1 PCS

6.2406.000

Mercury drop catcher

Silver wire (to capture mercury drops by amalgamation)



---

1 PCS

6.2711.030

Tray for filling the MME

For 884 Professional VA, 797 VA Computrace, 757 VA Computrace and 647, 663, 694, 747 VA Stands



1 PCS

6.2816.020

Syringe 10 mL with Luer connection

For various applications in IC and VA



---

1 PCS

6.2816.030

Needle with Luer connector

Needle for filling in mercury into the Multi-Mode Electrode and for sample addition in Karl Fischer titration.

