



Equipment: 807 Dosing Unit for pipetting system for Professional VA/CVS instruments

6.5339.620

Accessory kit for fully automated CVS pipetting system in combination with the 858 Professional Sample Processor. Contains one 807 Dosing Unit 5 mL and 50 mL each, clear glass bottles 2 L, pipetting tubing, rinsing station, safety shield and other accessories. Without Dosinos.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.5339.620

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.1456.210

Measuring vessel / 10-90 mL

Measuring vessel for automated VA systems. Dilution vessel for Advanced IC Sample Processor



1 PCS

6.1462.270

Robotic arm with transfer head, for 786 Swing Head, right-swinging

Robotic arm for pipetting system with 858 Professional Sample Processor VA and CVS



1 PCS

6.1562.120

Pipetting tubing / 10 mL with holder for sample processors

Pipetting tubing for liquid transfer.



1 PCS

6.1602.115

Bottle attachment / S 40 / green

For direct dosing from reagent bottles with S 40 thread (Merck...).





With light and kink protection.



With light and kink protection



With light and kink protection



1 PCS

6.1808.010

Coupling piece 2 X M6

2 M6 outer threads. Part of the screw connections for Exchange Units and 676 Sample Changer



1 PCS

6.1808.050

Tubing adapter olive / M8 outer

1 M8 outer thread and 1 tubing olive. E.g. for thermostat jacket of exchange units and stability measuring instruments.



1 PCS

6.1808.060

T connector / M6 inner

3 M6 inner threads.



1 PCS

6.1819.080

Cannula for gas feeding

For stability measuring instruments.



1 PCS

6.1823.010

Guide shaft to pipetting tubing

For assembly in the transfer robotic arm on the Robotic Sample Processors



2 PCS

6.1829.080

Aspiration tubing / M6 / 0.3 m



1 PCS

6.2622.040

Positioning plate

Positioning plate for pipetting system with 858 Professional Sample Processor VA and CVS



1 PCS

6.2622.030

Safety shield accessories

Accessories for pipetting system with 858 Professional Sample Processor VA and CVS. Contains stand and fastening material for the 6.2751.200 safety shield and the 6.2622.020 positioning plate.



2 PCS

6.2709.070

Guiding sleeve

For fixation of buret and dosing tips in titration heads and electrode holders



2 PCS

6.2744.200

Adapter UNF 10/32 outer / M6 inner

Connecting 1/16 in. PEEK capillaries and transfer tubing (M6)



1 PCS

6.2751.200

Safety shield

Safety shield for pipetting system with 858 Professional Sample Processor VA and CVS



1 PCS

6.2841.110

Rinsing station for pipetting equipment

Rinsing station for pipetting equipment



807 Dosing Unit with integrated data chip with 5 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread. FEP tubing connection, antidiffusion tip.



807 Dosing Unit with integrated data chip with 50 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread. FEP tubing connection, antidiffusion tip.



Optional accessories

Order no.	Description
2.800.0020	800 Dosino Drive with write/read hardware for intelligent dosing units. With permanently attached cable (length 0.65 m).

