



VA electrode shaft accessories with SPE electrode shaft for VA Computrace

6.5340.030

Accessory equipment for use with *screen-printed electrodes* (SPE). Contains electrode shaft for screen-printed electrodes, measuring vessel, stirrer, and additional accessories. Without electrodes.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.5340.030

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS 6.1204.200 Stirrer for VA Stand
Stirrer for VA Stand, complete, metal-free




1 PCS 6.1241.090 Electrode shaft for screen-printed electrodes (SPE)

Holder for *screen-printed electrodes* (SPE) with plug contacts for connecting the electrode cables for working, reference, and auxiliary electrodes. When used with screen-printed electrodes, the electrode shaft provides a complete sensor for voltammetry. Whereas the screen-printed electrodes must be replaced at regular intervals, the electrode shaft can continue to be used.



1 PCS 6.1244.020 Driving belt
For Rotating Disk Electrodes and stirrer in VA Stands (3 pieces).



1 PCS 6.1244.060 Silicone seal for SPE electrode shaft (10 pieces)

Seal made of silicone for SPE electrode shaft 6.1241.090 and measuring head (SPE) 6.01256.020, 10 pieces.



1 PCS

6.1415.150

Titration vessel / 5-70 mL



1 PCS

6.1446.000

SGJ stopper / B-14/(15)



4 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.1808.000

Connection bushing M6



1 PCS

6.1819.000

Gas inlet tube

In lid of measuring vessel in VA Stand



3 PCS

6.2112.130

Adapter plug contact 5.5 mm / screw contact, gold-plated

For VA electrodes

