



Consumable Kit Rancimat

6.5706.020

Assembly of important expendable items for the Rancimat.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.5706.020

Qt.	Order no.	Description
1 PCS	6.1428.107	Polystyrene measuring vessel for stability measurements

Conductivity measuring vessel, made of polystyrene, for stability measuring instruments. Set of 50 pieces.



1 PCS

6.1429.040

Reaction vessel short for stability measurements

Reaction vessel for stability measurements with Rancimat and PVC Thermomat. Set of 117 pieces.




2 PCS

6.1454.040

O-Ring

O-ring for air tube for stability measurements. Set of 8 pieces.



8 PCS

6.1816.010

Silicone tubing 0.22 m

Connection tubing for stability measuring instruments.



1 PCS

6.2418.100

Air tube for oil measurements

Set of 117 pieces. In combination with 6.1429.040.



2 PCS

6.2724.010

Dust filter

Dust filter for Rancimats and Karl Fischer ovens.



1 PCS

6.2753.107

Reaction vessel cover for stability measurements

For 6.1429.040 or 6.1429.050 Reaction vessel. Set of 100 pieces.



Optional accessories

Order no.	Description	
6.2811.000	Molecular sieve Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.	
6.2811.010	Molecular sieve (US) Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments. This product is intended for all users who are US residents. The material is rated hazardous in the US in accordance with OSHA Hazard Communication Standard (29 CFR 1910.1200).	