



NIRS XDS Process Analyzer – SingleFiber 4 Channels

A629280220C

The NIRS XDS Process Analyzer – SingleFiber provides the next generation of process analyzers for real-time analysis in the pharmaceutical and chemical industries. Non-destructive, accurate measurements of clear liquids are performed directly in the process line or in a reaction vessel. All 4 channels can be configured independently from each other in the software. Single fiber optics allows the 4 measuring points to be up to 100 meter away from the XDS NIR process analyzer.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
-----------	-------------

A649281190 NIRS XDS Process Contact Lamp 50W, 12 Volt- White lead

Replacement lamp with contact for NIRS Process instruments. 50W, 12 Volt- White lead.

