



NIRS Analyzer PRO – ContactReflection

A629281110

The NIRS Analyzer Pro is a single point process analysis system based on high resolution diode array technology. It provides non-destructive analysis of solids, granules, or powders through a glass window directly in the process line without bypass – a true inline system.


Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
6.7440.480	NIRS Process, I/O Communication box
	I/O Communication box for process instruments with 24 analog and 24 digital channels



A660044277

NIRS Analyzer PRO Lamp 12V/20W

Replacement lamp for the NIRS Analyzer Pro (12V/20W)

