



NIRS Analyzer PRO – DirectLight/NonContact

A629281130


The NIRS Analyzer Pro is a single point process analysis system based on high resolution diode array technology. It provides non-destructive analysis of solids like granules and powders. This system is typically mounted above a conveyor belt in order to analyze products passing by. Fast moving, uneven surfaces are able to be analyzed easily with this configuration.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.


Scope of delivery A629281130

Qt.	Order no.	Description
-----	-----------	-------------

I/O Communication box for process instruments with 24 analog and 24 digital channels



Optional accessories

Order no.	Description	
6.7440.480	NIRS Process, I/O Communication box I/O Communication box for process instruments with 24 analog and 24 digital channels	
A660044277	NIRS Analyzer PRO Lamp 12V/20W Replacement lamp for the NIRS Analyzer Pro (12V/20W)	