



TacticID-N Plus Handheld Raman Spectrometer

BWT-840000994


The TacticID[®]-N Plus is a field-ready handheld Raman Spectrometer specifically designed for non-contact forensic analysis of narcotics, pharmaceutical drugs, cutting agents, and their precursors by law enforcement personnel.

Featuring an intuitive workflow and touch screen, samples can be nondestructively analyzed through opaque and transparent packaging. The sample threat level is displayed prominently to the users, e.g. security staff, emergency services (such as law enforcement authorities), customs and border guards as well as bomb disposal squads, and hazmat teams, allowing them to act quickly with minimal sample contact. The TacticID-N Plus utilizes laboratory-proven Raman spectroscopy, which allows users to obtain identification of illicit substances in real time without ever compromising the integrity of the sample or the chain of evidence.

TacticID-N Plus users have access to regular library updates in order to continuously maintain up-to-date identification capabilities and stay ahead of emerging narcotics.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery BWT-840000994

Qt.	Order no.	Description
1 PCS	BWT-840000295	Vials for liquids, 15 mm in diameter Borosilicate glass vials (diameter 15 mm), package of 6, with pouch.
		
1 PCS	BWT-840000403	10.8 V rechargeable lithium ion battery Rechargeable lithium-ion battery: fits TacticID-N Plus, TacticID Mobile, and NanoRam-1064. TID-BAT
		
1 PCS	BWT-840000406	TacticID Vial-Holder Adapter TacticID Vial-Holder Adapter, compatible with 15 mm-diameter vials.
		
1 PCS	BWT-840000414	18 V Power Adapter Power adapter for 18 V DC, 1.67 A. Plugs for USA, EU, UK and AU included in the scope of delivery. TID-PWA
		

1 PCS

BWT-840000416 TacticID Point-and-Shoot adapter

TacticID Point-and-Shoot adapter.



1 PCS

BWT-840000417 TacticID Validation Cap

TacticID Validation Cap with integrated ASTM-standard polystyrene. For use with TacticID-N Plus and TacticID-GP Plus.



1 PCS

BWT-840000419 TacticID Stylus/pen combo

Stylus/pen combo for touch screens (package of 5).

TID-STY



1 PCS

BWT-840000422 TacticID Carrying Case

TacticID hard plastic rugged carrying case.



1 PCS

BWT-840000429 TacticID-N shaft protection cap (Black)

TacticID-N and TacticID-N Plus shaft protection cap (Black).

TID-CAP-N



1 PCS

BWT-840000858 USB cable

USB cable supporting micro USB port up to type A USB 2.0.



1 PCS

BWT-840000948 USB OTG Host Cable

USB OTG host cable for data transmission between a USB flash drive and a handheld system.



1 PCS

BWT-840000981 TID-SFW-PS TacticID Plus Software

USB flash drive containing TacticID Plus TID software for PC, TID installation package and TacticID Plus Product Manuals.



TactiID N Plus library of illicit materials includes drugs, TICs, TIMs, precursors, and general chemicals.



TAP-KIT-T is a SERS Analysis kit for use with the TactiID-N Plus for Raman signal amplification and fluorescence reduction of hard-to-measure samples. This kit includes a black pouch, Tweezer, SERS Analysis Adaptor, and 25 SERS substrates

SERS analysis kit for use with Raman handheld spectrometers from B&W Tek. For Raman signal amplification and fluorescence reduction of hard-to-measure samples. For use with 785 nm laser excitation. Includes: tweezers, TacPac adapter, 25 SERS strips.

TAP-KIT-T




Pouch for handheld Raman adapters compatible with NanoRam-1064, TactiID-N Plus, and TactiID Mobile.

NR2-PCH-ADPT



Optional accessories

Order no.	Description	
BWT-840000294	Battery Charging Cradle	
	Lithium-ion battery charging cradle.(for charging NR-BAT and TID-BAT only).	
BWT-840000329	Pouch for 15 mm Vials	
	Pouch for 15 mm vials.	
BWT-840000353	Shaft Protection Cap	
	Handheld Raman shaft protective cap, package of 10.	
BWT-840000361	Handheld Laser Power Meter	
	Laser Power Meter, Pen Type, Calibrated. Range: 5 micro-Watts to 1 Watt. Used in OQ/PQ and annual certification laser power measurement.	

BWT-840000912 O-ring replacement kit for handheld Raman instruments

Handheld Raman O-ring replacement kit –(for shaft only).

Kit includes: O-Ring#1 ID 9.00 mm x15 and O-Ring#2 ID 10.00 mm x15 (TacticID-N Plus, TacticID Mobile, and NanoRam-1064).

TID-ORK



BWT-840000491 TacticID Right-Angle Adapter

TacticID Right-Angle Adapter.



BWT-840000498 Disposable Battery Cartridge

TacticID battery cartridge for use with six type 123A disposable batteries (e.g. Surefire). Batteries not included.



BWT-840000524 TacPac™ Adapter


SERS analysis adapter for use with TacPac™-P SERS substrates.



TactiD (785 nm) Pouch for Adapters



Car charger for TactiD-N Plus.



Raman Spectral Library for TactiD Plus with more than 1,180 spectra of common chemicals, narcotics, solvents and plastics. Includes NFPA704 and GHS safety instructions.



SERS analysis sample preparation kit for use with Raman handheld spectrometers from B&W Tek. For Raman signal amplification and fluorescence reduction of hard-to-measure samples. For use with 785 nm laser excitation. Disposable plastic single-use laboratory spatula, pipettes and PP graduated microcentrifuge tubes (2 ml) used to prepare samples for testing with TacPac™ SERS. Package of 25 SERS substrates included.

TAP-SAM25



BWT-840001125 TacPac sample preparation kit (package of 50)

SERS analysis sample preparation kit for use with Raman handheld spectrometers from B&W Tek. For Raman signal amplification and fluorescence reduction of hard-to-measure samples. For use with 785 nm laser excitation. Disposable plastic single-use laboratory spatula, pipettes and PP graduated microcentrifuge tubes (2 ml) used to prepare samples for testing with TacPac™ SERS. Package of 50 SERS substrates included.



TAP-SAM50

607506480 SERS Discovery Kit

Introductory kit for SERS analysis containing gold and silver P-SERS strips and gold and silver colloids.

