

Cell for Screen-Printed Electrodes in PEEK

DRP-CELL-PEEK


Cell in PEEK for batch analysis which allows optional stirring and it's recommended for deaerated samples and for standard additions mode. Suitable for working with all catalogue Screen-Printed Electrodes but recommended to be used with work in solution format SPEs.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description	
DRP-CAC	Cable connector for Screen-Printed Electrodes	
	Connects DropSens Screen-Printed Electrodes to any kind of potentiostat. Recommended when working in batch or flow analysis.	
DRP-CAST	µStat Cable connector for Screen-Printed Electrodes	
	Connects any DropSens equipment with DropSens Screen-Printed Electrodes	
DRP-CELLHOLDER	Cell Holder	
	Accessory to be used in combination with CELL and CELL-PEEK when working in batch analysis for set up improvement.	