



Flow Cell in Teflon - Only Working Electrode

DRP-CFLWCL-WE-TEF

Teflon wall-jet Flow-Cell for FIA. Its o-ring entraps only the working electrode. Suitable to be used with standard Screen-Printed Electrodes with the electrochemical cell in the center of the strip. Closing system with powerful magnets. Fittings included.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
-----------	-------------


Polyaniline modified Screen-Printed Carbon Electrode designed for the development of (bio) sensors with an enhanced electrochemical active area.



DRP-CAC

Cable connector for Screen-Printed Electrodes

Connects DropSens Screen-Printed Electrodes to any kind of potentiostat. Recommended when working in batch or flow analysis.



DRP-710-U75

Screen-Printed Prussian Blue/Carbon Electrode

Screen-Printed Prussian Blue/Carbon Electrode. Ideal for the determination of hydrogen peroxide at a low detection potential.



DRP-CAST

μ Stat Cable connector for Screen-Printed Electrodes

Connects any DropSens equipment with DropSens Screen-Printed Electrodes



DRP-550-U75

Screen-Printed Platinum Electrode (Aux.:Pt; Ref.:Ag)

Screen-Printed Platinum Electrode (Aux.:Pt; Ref.:Ag). Suitable for working with microvolumes, for decentralized assays or to develop specific sensors.

