




Flow Cell in Teflon for Interdigitated Electrodes


DRP-FLWCL-IDE-TEF

Teflon wall-jet Flow-Cell for FIA. Suitable to be used with IDE electrodes ref. G-IDEAU10, G-IDEAU5, G-IDEPT10, G-IDEPT5, G-IDECONAU10 and G-IDECONPT10. Closing system with powerful magnets. Fittings included.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description
DRP-CACIDE	Cable connector for Interdigitated Electrodes/Glass Substrate
	Connects DropSens Interdigitated Electrodes over Glass Substrate with any kind of potentiostat



DRP-G-IDEAU10-U20

Interdigitated Gold Electrode / 10 microns lines and gaps / Glass Substrate

Interdigitated Gold Electrode / 10 microns lines and gaps / Glass Substrate



DRP-G-IDECONAU10-U20

Interdigitated Gold Concentric Electrode / 10 microns lines and gaps / Glass Substrate

Interdigitated Gold Concentric Electrode / 10 microns lines and gaps / Glass Substrate



DRP-G-IDEPT5-U20

Interdigitated Platinum Electrode / 5 microns lines and gaps / Glass Substrate

Interdigitated Platinum Electrode / 5 microns lines and gaps / Glass Substrate

