



NIRS XDS Interactance OptiProbe Module

2.921.0510

Robust measurement system for the monitoring of reactions in technical centers and in pilot plants

The NIRS XDS Interactance OptiProbe Analyzer enables the simple and reliable monitoring of chemical reactions in technical centers and pilot plants. The method development and the scale-up of production processes are areas of application in which the NIRS XDS Interactance OptiProbe Analyzer supplies precise results regarding the identity and quality of a wide variety of samples. Whereas solids, strongly scattering liquids and slurries are measured with the reflection probe, the immersion probe is used to analyze aqueous products, clear liquids and solvents. The probe is connected to the analyzer through an optical fiber to enable secure measurements even under rough and /or dangerous conditions in the processing environment.

納品範囲 2.921.0510

Qt.

Order no.

説明

Robust measurement system for the monitoring of reactions in technical centers and in pilot plants

The NIRS XDS Interactance OptiProbe Analyzer enables the simple and reliable monitoring of chemical reactions in technical centers and pilot plants. The method development and the scale-up of production processes are areas of application in which the NIRS XDS Interactance OptiProbe Analyzer supplies precise results regarding the identity and quality of a wide variety of samples. Whereas solids, strongly scattering liquids and slurries are measured with the reflection probe, the immersion probe is used to analyze aqueous products, clear liquids and solvents. The probe is connected to the analyzer through an optical fiber to enable secure measurements even under rough and/or dangerous conditions in the processing environment.



オプションアクセサリ

Order no.	説明
6.7430.080	XDS Interchangeable Probe Set 交換用プローブセットには、反射プローブ、半透過プローブ、様々な長さ(2cm、4cm、5cm、10cm)に設定できる固定用ねじを含むスペーサーが含まれています。更に、プローブの保護カバー、六角レンチ、パッキングリングが付属しています。 このセットはXDS SmartProbe AnalyzerおよびXDS Interactance Optiprobe Analyzerの付属品です。 The image shows the components of the XDS Interchangeable Probe Set. It includes two long, thin, cylindrical metal probes, one of which is shown with a protective cap. There are also four small, cylindrical metal spacers of different lengths, a hex key (Allen key), and a small metal ring (packing ring).