





Flow Cell in Teflon for Interdigitated Electrodes

DRP-FLWCL-IDE-TEF

オプションアクセサリ

Order no.	説明
DRP-CACIDE	Cable connector for Interdigitated Electrodes/Glass Substrate



DRP-G-
IDEAU10-U20

Interdigitated Gold Electrode / 10 microns lines and gaps / Glass
Substrate



DRP-G-
IDECONAU10-
U20

Interdigitated Gold Concentric Electrode / 10 microns lines and
gaps / Glass Substrate



DRP-G-IDEPT5-
U20

Interdigitated Platinum Electrode / 5 microns lines and gaps /
Glass Substrate

