





Electrochemical Quartz Crystal Microbalance module

EQCM.S

納品範囲 EQCM.S

Qt.	Order no.	説明
2 PCS	EQCM.AU	EQCM 6 MHz Au/TiO ₂ crystal working electrode



1 PCS

EQCM.CE

EQCM counter electrode



1 PCS

EQCM.REF.EL.S

Reference Electrode Ag/AgCl 3M KCl



オプションアクセサリ

Order no.	説明
-----------	----

EQCM.AU	EQCM 6 MHz Au/TiO2 crystal working electrode
---------	--



EQCM.PT	EQCM 6 MHz Pt-coated crystal working electrode
---------	--

