



OMNIS Solvent Module

2.1009.0010

OMNIS Solvent Module for reliable addition of solvent and aspirating off titration vessel content without contact with the reagents. New, improved version resistant against chloroform, chlorobenzene, xylene, and glacial acetic acid (requires OMNIS 4.4.1 or higher). Manual control or actuation via OMNIS Software for maximum user-friendliness. Including accessories, solvent, and waste bottle.

Phạm vi giao hàng 2.1009.0010

Qt.	Order no.	Description
1 PCS	1.1009.0010	OMNIS Solvent Module

OMNIS Solvent Module for reliable addition of solvent and aspirating off titration vessel content without contact with the reagents. Manual operation or control via OMNIS Software for maximum user-friendliness. Including accessories, solvent and waste bottle.



1 PCS

6.01543.000

M8 aspiration tip

For aspirating off reagent in connection with the OMNIS Solvent Module.



1 PCS

6.01600.200

Siphon Breaker

Prevents continued running of liquid out of the reagent bottle into the titration cell.




2 PCS

6.01804.210

PVC tubing 4/6 0.35 m

For tubing connections to the OMNIS Solvent Module or the Solvent Pump.



1 PCS

6.02801.000

Cotton for adsorber tube

Wad of cotton (4 pieces), suitable for OMNIS 601406010 adsorber tube.



1 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.1446.090

Stopper / B-14/15 / M8



1 PCS

6.1543.110

Dosing tip / M8 thread

Transfer tip with M8 thread



1 PCS

6.1602.105

Bottle attachment / GL 45 / green

For direct dosing from reagent bottles with GL 45 thread.



1 PCS

6.1608.023

Amber glass bottle / 1000 mL / GL 45

For Exchange Units. Bottle for auxiliary solutions.



1 PCS

6.1608.030

Round glass bottle / 1000 mL / GL 45



1 PCS

6.1623.000

Bottle cap overflow protection

Bottle cap overflow protection for 1 waste bottle, prevents the waste bottle from overflowing when an aspiration pump is used



2 PCS

6.1805.200

PTFE tubing / M8 / 0.5 m



1 PCS

6.1808.050

Tubing adapter olive / M8 outer

1 M8 outer thread and 1 tubing olive. E.g. for thermostat jacket of exchange units and stability measuring instruments.



1 PCS

6.1819.050

Tube

For M8 threaded holes of 6.1602.10X Bottle attachments with 6.1608.030 Aspiration bottle. For 703 and 803 Ti Stand, 681 Pump Unit and 679 Rancimat



1 PCS

6.2739.000

Wrench

For tightening connectors



1 PCS

6.2811.000

Molecular sieve

Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.



1 PCS

6.01807.000

Adsorption cartridge

Adsorption cartridge for OMNIS Solvent Module or Solvent Pump, for filling with 6.2811.000 molecular sieve.



1 PCS

6.02004.030

Bottle holder

Bottle holder for 2.1009.0010 OMNIS Solvent Module.



1 PCS

6.01805.130

3/250 mm PTFE tubing with screw nipple

Aspiration tubing with screw nipple to connect the reagent bottle to the Siphon Beaker.



Phụ kiện tùy chọn

Order no.	Description	
6.2811.000	Molecular sieve Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.	
6.2811.010	Molecular sieve (US) Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments. This product is intended for all users who are US residents. The material is rated hazardous in the US in accordance with OSHA Hazard Communication Standard (29 CFR 1910.1200).	
6.1446.220	Reagent Replacement Solvent Module stopper Stopper for reagent replacement in a coulometric titration cell with a Solvent Module.	