



OMNIS Coulometer without stirrer

2.1018.0010

OMNIS Coulometer without magnetic stirrer for stand-alone operation or as the heart of an OMNIS titration system. The coulometer can be enhanced with additional titration, dosing or coulometer modules.

Coulometry is the ideal method for water content determination in the trace range (10 µg to 10 mg absolute water content) in liquids, solids and gases. In addition, coulometric Karl Fischer titration is an absolute method making titer determination unnecessary.

Phạm vi giao hàng 2.1018.0010

Qt.	Order no.	Description
1 PCS	1.1018.0010	OMNIS Coulometer without stirrer
		OMNIS Coulometer without magnetic stirrer



Function license for OMNIS Coulometer, includes

- Liquid Handling (LQH)
 - Coulometric Karl Fischer titration (KFC)
 - Coulometric determination of the bromine index according to ASTM D1492 (BRC)
 - MEAS U/T/pH/CONC
-



Phụ kiện tùy chọn

Order no.	Description	
6.05406.000	KFC equipment for OMNIS Coulometer with generator electrode with diaphragm Complete accessory set for coulometric Karl Fischer titration.	
6.05406.010	KFC equipment for OMNIS Coulometer with generator electrode without diaphragm Complete accessory set for coulometric Karl Fischer titration.	
2.1009.0010	OMNIS Solvent Module OMNIS Solvent Module for reliable addition of solvent and aspirating off titration vessel content without contact with the reagents. New, improved version resistant against chloroform, chlorobenzene, xylene, and glacial acetic acid (requires OMNIS 4.4.1 or higher). Manual control or actuation via OMNIS Software for maximum user-friendliness. Including accessories, solvent, and waste bottle.	
6.1446.220	Reagent Replacement Solvent Module stopper Stopper for reagent replacement in a coulometric titration cell with a Solvent Module.	

6.02102.010 Cable MDL PL/SO 0.5 m

Connecting cable for OMNIS instruments with MDL connector (0.5 m)



6.02102.020 Cable MDL PL/SO 1 m

Connecting cable for OMNIS instruments with MDL connector (1 m)



6.02102.030 Cable MDL PL/SO 2 m

Connecting cable for OMNIS instruments with MDL connector (2 m)



6.02102.040 Cable MDL PL/SO 5 m

Connecting cable for OMNIS instruments with MDL connector (5 m)



6.02102.050 Cable MDL PL/SO 10 m

Connecting cable for OMNIS instruments with MDL connector (10 m)



6.06001.100 Function license DET

"Dynamic equivalence point titration" function license for one OMNIS titration system.



6.06001.110 Function license MET

"Monotonic equivalence point titration" function license for one OMNIS titration system.



6.06001.160 Function license KFT with conditioning

Function license for volumetric Karl Fischer titration with conditioning for an OMNIS titration system.



6.06001.320 Function license for 5 parallel titrations

Function license for a maximum of five parallel titrations using one OMNIS Titrator with additional Titration Modules or Dosing Modules.

