

733 IC Separation Center

2.733.0120

733 IC Separation Center with two injectors

Phạm vi giao hàng 2.733.0120

Qt.	Order no.	Description	
1 PCS	6.1816.000	Silicone tubing 6 mm, 0.5 m Drainage tube for IC Separation Center chamber	
1 PCS	6.1831.010	PEEK capillary 0.25 mm i.d. / 3 m For all IC components.	
2 PCS	6.2027.030	Column holder 8.5 mm For mounting columns in column chamber of IC instruments.	

2 PCS

6.2027.040

Column holder 11.3 mm

For mounting columns in column chamber of IC instruments.



2 PCS

6.2744.010

Pressure screw 5x

With UNF 10/32 connection. For the connection of PEEK capillaries



2 PCS

6.2744.020

Luer/UNF coupling

For IC instruments



4 PCS

6.2744.040

2 x UNF 10/32 coupling

For connecting 1/16 in. capillaries. For IC instruments



2 PCS


6.2816.020

Syringe 10 mL with Luer connection

For various applications in IC and VA



Phụ kiện tùy chọn

Order no.	Description	
6.1825.210	PEEK sample loop 20 μ L For injection valve, with 2 PEEK pressure screws	
6.2057.010	Holder for dialysis cell For the dialysis cell in IC column chambers	
6.5904.030	Rotor, 6/2 injector Rotor for the 6/2 injector in IC systems.	
6.5904.050	Injector with 0.25 μ L internal loop Flat plate injector with 0.25 μ L internal loop	

6.5904.060

Rotor (internal loop, 0.25 μ L)

Rotor for 4-port flat plate injector (6.5904.050), 0.25 μ L internal loop



6.5904.070

Rotor (internal loop, 1.0 μ L)

Rotor for 4-port flat plate injector (6.9956.001), 1.0 μ L internal loop

