




Sample changer KFT 760

2.760.0020

Complete system for automatic volumetric Karl Fischer titration.

Phạm vi giao hàng 2.760.0020

Qt.	Order no.	Description
1 PCS	6.2752.010	Plug cover
		Used with the 717, 730, 760 Sample Changers, IC Sample Processors, 774 Oven Sample Processor




1 PCS

1.741.0010

741 Magnetic Stirrer

Magnetic stirrer for sample changers.

The magnetic stirrer is located beneath the sample rack.



1 PCS

6.0340.000

Double Pt-wire electrode for sample changer (fixed cable)

Indicator electrode (diameter = 5.3 mm, fixed cable 2 m) used for automated Karl Fischer titration.



24 PCS

6.1432.210

Sample beaker, clear glass, 75 mL, 1 piece

Sample beaker, clear glass, 75 mL, suitable for sample racks: 6.2041.340, 6.2041.350, 6.2041.800 and 6.2041.900, 1 piece



1 PCS

6.1458.030

KF titration head / 4 x M10 connectors

With 4 M10 screw nipples, including stoppers.



2 PCS 6.1805.110 FEP tubing / M6 / 80 cm
With light and kink protection



1 PCS 6.1815.010 Spiral band 0.5 m
For holding together different cables or tubing.



1 PCS 6.1821.000 Aspiration tube
For KF Sampler



24 PCS 6.1903.010 Stirring bar / 12 mm
Stirring Bar with magnetic core, PTFE coated, length 12 mm



24 PCS

6.2037.040

Foil holder

Foil holder for 75 mL sample beaker for Karl Fischer for KF Sample Changer



1 PCS

6.2041.340

Sample rack, 24 x 75 mL, (PVC)

Sample rack for 810 Sample Processor, 814 USB Sample Processor, 855 Robotic Titrosampler "Light", suitable for 24 sample beakers. The following sample beakers can be used: 6.1432.210, 6.1459.400.

Plastic: Polyvinyl chloride (PVC)



1 PCS

6.2820.000

Al foil

1000 circular sheets of 80 mm diameter. Thickness: 0.010 mm. Used together with 6.2037.030 Foil holder for 6.1432.210 Sample beaker. For 674 KF Sample Changer

