



786 Swing Head with transfer head left

2.786.0010

Additional equipment for Robotic Sample Processors XL for transferring / pipetting from small sample vessels into larger titration vessels on the rack or into an external titration cell; this can be located to the left of the work station.

Phạm vi giao hàng 2.786.0010

Qt.	Order no.	Description
1 PCS	1.786.0010	786 Swing Head with transfer head left
		Additional equipment for Robotic Sample Processors XL for transferring / pipetting from small sample vessels into larger titration vessels on the rack or into an external titration cell; this can be located to the left of the work station.



1 PCS

6.1462.030

Robotic arm with transfer head for 786 Swing Head,
left swinging

Transfer robotic arm for liquid handling purposes



1 PCS

6.1562.100

Pipetting tubing / 10 mL

Pipetting tubing for fully automatic transfer of sample solutions.



1 PCS

6.1823.010

Guide shaft to pipetting tubing

For assembly in the transfer robotic arm on the Robotic Sample Processors



1 PCS

6.2621.120

Hexagon key 1.5 mm

