



805 Dosimat

2.805.0010

Dosing device for the Titrando and Sample Processors with read/write hardware for intelligent Exchange Units. With permanently attached cable. Without Exchange Unit.

Phạm vi giao hàng 2.805.0010

Qt.	Order no.	Description
1 PCS	1.805.0010	805 Dosimat
		Dosing device for the Titrando and Sample Processors with read/write hardware for intelligent Exchange Units. With permanently attached cable. Without Exchange Unit.



1 PCS

6.2739.010

Wrench

For exchange units.



Phụ kiện tùy chọn

Order no.	Description	
6.2151.010	Cable MSB plug / socket Extension cable for MSB connections and for Touch Control - Titrandos, Mini-DIN socket - mini-DIN plug	
6.3026.110	Exchange Unit 1 mL Exchange Unit with integrated data chip with 1 mL glass cylinder and light protection. PCTFE/PTFE flat stopcock, FEP tubing connection, antidiffusion buret tip and standard amber glass reagent bottle.	
6.3026.150	Exchange Unit 5 mL Exchange Unit with integrated data chip with 5 mL glass cylinder and light protection. PCTFE/PTFE flat stopcock, FEP tubing connection, antidiffusion buret tip and standard amber glass reagent bottle.	
6.3026.210	Exchange Unit 10 mL Exchange unit with integrated data chip with 10 mL glass cylinder and light protection. PCTFE/PTFE flat stopcock, FEP tubing connection, antidiffusion buret tip and standard amber glass reagent bottle.	

6.3026.220

Exchange Unit 20 mL

Exchange unit with integrated data chip with 20 mL glass cylinder and light protection. PCTFE/PTFE flat stopcock, FEP tubing connection, anti-diffusion buret tip and standard amber glass reagent bottle.



6.3026.250

Exchange unit 50 mL

Exchange unit with integrated data chip with 50 mL glass cylinder and light protection. PCTFE/PTFE flat stopcock, FEP tubing connection, antidiffusion tip and standard amber glass reagent bottle.

