

813 VA Compact Autosampler


2.813.0020

Autosampler for the basic automation of the VA systems, with integrated peristaltic pump for sample transfer. The sample changer can be used in combination with the 746 VA Trace Analyzer or the 797 VA Computrace. Capacity: 18 samples.

Phạm vi giao hàng 2.813.0020

Qt.	Order no.	Description
1 PCS	6.2751.070	Safety shield

Safety shield for the 813 Compact IC Autosampler.




1 PCS	6.2752.010	Plug cover
-------	------------	------------

Used with the 717, 730, 760 Sample Changers, IC Sample Processors, 774 Oven Sample Processor



1 PCS	6.1822.410	PEEK connection capillary
-------	------------	---------------------------

Connection capillary for sample transfer with sample changers, 1 x M6, inner diameter 0.5 mm, for trace analysis




2 PCS

6.1826.020

PVC pump tubing (blue/blue)

For online IC instruments and automation in voltammetry.



1 PCS

6.1831.050

PEEK capillary 0.5 mm i.d., 40 cm

For inline dialysis and for sample changers. For IC Dialysis Unit, IC Sample Processor, Compact IC Autosampler, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit, Compact VA Autosampler



1 PCS

6.2621.060

5/16 in. wrench

For 5/16 in. screws. For IC instruments



2 PCS

6.2621.090

1/2 in. wrench

For 1/2 in. screws. For IC Sample Processors





200 pieces. For IC Sample Processors and VA Autosamplers.

The IC sample vessels are made of polypropylene (PP) plastic for health-care applications. Production and packaging both take place under sterile and clean room conditions (Clean Room Class 7, in accordance with DIN ISO 14644). This guarantees minimal degrees of contamination. The sample vessels are subjected to random checks for leachable anions and cations.



With UNF 10/32 connection. For the connection of PEEK capillaries



Connection between pressure screw and pump tubing. 4 pieces. For IC instruments with peristaltic pump



1 PCS

6.2833.000

Needle/tube holder 1/8 in.

For IC/VA Sample Processors. For 1/8 in. tubes



Phụ kiện tùy chọn

Order no.	Description	
6.2141.170	Cable 797 VA Computrace to 813 Compact Sample Changer to 731 Relay Box	
2.772.0110	772 Pump Unit Peristaltic pump for direct connection to Sample Processors.	
2.823.0010	823 Membrane Pump Unit Membrane pump for direct connection to Sample Processors.	