



843 Pump Station (membrane) - rinse/aspirate for Compact Sample Changer

2.843.0030

The 843 Pump Station (membrane) has two built-in membrane pumps. These can be controlled directly from the 869 Compact Sample Changer via remote signals. The Rinse / Aspirate version is provided with all the accessories needed for automatically emptying the titration beaker and rinsing the titration equipment. The accessories is defined for usage with the 869 Compact Sample Changer.

Phạm vi giao hàng 2.843.0030

Qt.	Order no.	Description
1 PCS	1.843.0020	843 Pump Station (membrane)
The 843 Pump Station (membrane) has two built-in membrane pumps. These can be controlled directly via the interface using remote signals or manually by pressing a button.		



1 PCS

6.1543.120

Aspiration tip / M8 thread

Aspiration tip if 6.1805.200 siphoning tubing is used



2 PCS

6.1621.000

Container 10 L

As rinsing or waste container in automated systems.



2 PCS

6.1621.100

Cap for 10-L container

Cap for 10-L container. For easy and safe transportation of filled containers.



4 PCS

6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m



1 PCS

6.1805.200

PTFE tubing / M8 / 0.5 m



3 PCS

6.1805.420

FEP tubing / M6 / 48 cm

With light and kink protection



1 PCS

6.1808.240

Distributor Rinsing / Aspiration for 869 Compact
Sample Changer

Distributor with 3 x M6, 1 x M8 and 2 x 1/8 connections for the use of rinsing and / or aspiration equipment on the 869 Compact Sample Changer.



2 PCS

6.1812.000

PTFE tubing 4 mm / 6 mm / 4 m



2 PCS

6.1828.000

PVDF connection nipple

For 6.1621.000 canister.



3 PCS

6.2740.020

Spray nozzle

For the fine-spraying of the rinsing solution.



Phụ kiện tùy chọn

Order no.	Description	
6.2141.260	Remote cable Compact Sample Changer - Dosimat plus and Titrino plus Cable for remote connection between Compact Sample Changer and Dosimat plus and further devices (e.g. Titrino plus or Pump Station).	
6.2141.300	Remote cable 843 Pump Station to Sample Processor Remote cable for the direct connection of the 843 Pump Station to the pump connectors of the Sample Processors	
6.2740.030	Rinsing nozzle The rinsing solution is applied as a fine jet for the optimal removal of layers on the electrodes and titration accessories. For 6.5613.000 ECO titration equipment.	