



889 IC Sample Center

2.889.0010

The 889 IC Sample Center is a robust Autosampler for high sample throughput, optimized for the challenges of the modern analytic laboratory. It works in accordance with the x-y-z principle and a syringe control with high resolution for precise sample injection. A double needle system enables the penetration of vessel caps and septa. A PEEK injection valve completes the automation system of the 889 IC Sample Center. It is the appropriate automation solution when you have only a very small amount of sample.

Phạm vi giao hàng 2.889.0010

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

1.889.0010

889 IC Sample Center

The 889 IC Sample Center is a robust Autosampler for high sample throughput, optimized for the challenges of the modern analytic laboratory. It works in accordance with the x-y-z principle and a syringe control with high resolution for precise sample injection. A double needle system enables the penetration of vessel caps and septa. A PEEK injection valve completes the automation system of the 889 IC Sample Center. It is the appropriate automation solution when you have only a very small amount of sample.



1 PCS

6.1608.110

Wash bottle 250 mL to IC Sample Center

Bottle for washing solution in the IC Sample Center.



2 PCS

6.1807.010

Y-connector for tubing 6-9 mm i.d.

Connector for waste tubing



2 PCS

6.1816.070

Silicone tubing 8/12 mm, 1 m

For the IC Sample Center



1 PCS

6.1831.010

PEEK capillary 0.25 mm i.d. / 3 m

For all IC components.



2 PCS

6.2041.200

Sample rack, 48 x 11 mL, (Al)

Sample rack for 889 IC Sample Center, suitable for 48 sample vessels. The following sample vessels can be used: 6.2743.100, 6.2743.107, 6.2743.110, 6.2743.117.

Material: Aluminum (Al)



1 PCS

6.2151.020

Cable USB A - USB B 1.8 m

USB connecting cable



1 PCS

6.2743.107

PP Vials 0.7 mL to IC Sample Center, 100 pieces

100 sample vials.



1 PCS

6.2743.127

Septum to 6.2743.1XX Si/PTFE

100 pieces



1 PCS

6.2744.010

Pressure screw 5x

With UNF 10/32 connection. For the connection of PEEK capillaries



1 PCS



6.2816.200

Air needle 80 mm, black

To IC Sample Center




Phụ kiện tùy chọn

Order no.	Description	
6.1562.170	Tefzel transfer tubing 1000 µL To IC Sample Center.	
6.1829.060	Capillary set for syringe unit For the IC Sample Center	
6.2743.100	PP Vials 0.7 mL to IC Sample Center, 1000 pieces 1000 sample vials.	
6.2743.110	PP Vials 0.3 mL to IC Sample Center, 1000 pieces 1000 sample vials.	

6.2743.117 PP Vials 0.3 mL to IC Sample Center, 100 pieces
100 sample vials.



6.2743.120 Septum to 6.2743.1XX Si/PTFE
1000 pieces



6.2816.210 Needle with capillary, 15 μ L
For the IC Sample Center



6.5333.000 IQ/OQ Kit for IC
The IQ/OQ Kit contains all parts and standard solutions required for IQ /OQ in ion chromatography.




6.7303.100 Glass vial, 1.5 mL, 100 pieces

For use with the following sample changers: Liquid Autosampler CIC (TEI), MMS 5000 Autosampler (AJ) and 889 Sample Center.



6.7303.110 Septum for 67303100, 100 pieces

For use with the 67303100 Glass vial on the following sample changers: Liquid Autosampler CIC (TEI), MMS 5000 Autosampler (AJ), and 889 Sample Center.



6.2041.210 Sample rack, 12 x 6 mL, (Al)


Sample rack for 889 IC Sample Center, suitable for 12 sample vessels. The following sample vessels can be used: 6.2419.000.

Material: Aluminum (Al)



6.1448.050 Aluminum septum caps / 1000 pieces

Including Al sealing.



6.2419.000 Sample glass, 6 mL, 1000 pieces

Sample glass for up to 6 mL sample, suitable for 6.1448.050 septum, 1,000 pieces



6.2621.110 Crimping tongs

Crimping tongs for sealing 6 mL sample vials (6.2419.000, 6.2419.007).



6.2420.007 Sample glass with thread, 8 mL, 100 pieces

Sample glass with thread, 8 mL, suitable for 6.1448.067 screw cap, 100 pieces



6.1448.067 Screw cap with septum, 100 pcs.

Screw cap with septum, suitable for sample glass 8 mL 6.2420.007, 100 pieces

