



NIRS XDS RapidContent Module

2.921.0110

Rapid, nondestructive analyses of solid and liquid substances

The NIRS XDS RapidContent Analyzer enables rapid, nondestructive analyses of solid and liquid formulations and substances. The sample is centered with the aid of an iris, thus enhancing the reproducibility of the analysis. A gold diffuse reflector enables the measurement of liquid samples. The Solids Module of the NIRS XDS RapidContent Analyzer extends the area of application to include practically every form of solid samples, whether it be fine powders, coarse granulates, pellets or flakes.

Phạm vi giao hàng 2.921.0110

Qt.	Order no.	Description
-----	-----------	-------------

Rapid, nondestructive analyses of solid and liquid substances

The NIRS XDS RapidContent Analyzer enables rapid, nondestructive analyses of solid and liquid formulations and substances. The sample is centered with the aid of an iris, thus enhancing the reproducibility of the analysis. A gold diffuse reflector enables the measurement of liquid samples. The Solids Module of the NIRS XDS RapidContent Analyzer extends the area of application to include practically every form of solid samples, whether it be fine powders, coarse granulates, pellets or flakes.



Iris for focusing mini sample vessels, reflection standard, various vials or a transflection vessel in the NIRS XDS RCA Analyzer, the NIRS XDS MasterLab Analyzer and the NIRS XDS MultiVial Analyzer



Phụ kiện tùy chọn

Order no.	Description	
6.7402.200	Stopper for 6.7402.110 123 stoppers for 15 mm Reflection disposable vials.	
6.7402.210	Stopper for 6.7402.120 225 stoppers for 19 mm Reflection disposable vials.	
6.7402.220	Stopper for 6.7402.130 172 stoppers for 22 mm Reflection disposable vials.	
6.7402.230	Stopper for 6.7402.140 216 stoppers for 28 mm Reflection disposable vials.	

6.7402.110

Disposable vials, 15 mm, reflection

123 lockable disposable glass vials with a diameter of 15 mm for analyses of solids in reflection. Suitable for NIRS RapidContent (Solid) Analyzers of the XDS, DS2500 and OMNIS product families.



6.7402.120

Disposable vials, 19 mm, reflection

225 lockable disposable glass vials with a diameter of 19 mm for analyses of solids in reflection. Suitable for NIRS RapidContent (Solid) Analyzers of the XDS, DS2500 and OMNIS product families.



6.7402.130

Disposable vials, 22 mm, reflection

172 lockable disposable glass vials with a diameter of 22 mm for analyses of solids in reflection. Suitable for NIRS RapidContent (Solid) Analyzers of the XDS, DS2500 and OMNIS product families.



6.7402.140

Disposable vials, 28 mm, reflection

216 lockable disposable glass vials with a diameter of 28 mm for analyses of solids in reflection. Suitable for the following analyzers:

- NIRS DS2500 Analyzer
- NIRS XDS RapidContent (Solid) Analyzer
- NIRS XDS MultiVial Analyzer
- NIRS XDS MasterLab Analyzer

