



943 Professional Reactor Vario

2.943.0110

The 943 Professional Reactor Vario is a heatable reactor, especially developed for pre- and post-column derivatization. It captivates by its robustness, the fast heating rate and the possibility of carrying out reactions at up to 150 °C.

Phạm vi giao hàng 2.943.0110

Qt.	Order no.	Description
1 PCS	1.943.0010	943 Professional Thermostat / Reactor Vario



1 PCS

6.1015.100

Metrosep BP 1 Guard/2.0

The Metrosep BP 1 Guard/2.0 is used to generate a sufficiently high working pressure in the flow path of post-column reagents. They are used in combination with a high-pressure pump for conveying the post-column reagent in the Professional Reactor.



1 PCS

6.1803.150

PTFE Capillary i.d. 0.25 mm / 1 m



2 PCS

6.1831.100

PEEK capillary 0.25 mm i.d., 1 m



1 PCS

6.2151.020

Cable USB A - USB B 1.8 m

USB connecting cable



1 PCS

6.2621.100

3 mm hex key for IC Sample Processors



1 PCS

6.2621.140

Hexagon key 2.5 mm



1 PCS

6.2744.010

Pressure screw 5x

With UNF 10/32 connection. For the connection of PEEK capillaries



1 PCS

6.2744.014

Pressure screw 2x

With UNF 10/32 connection. For the connection of PEEK capillaries



1 PCS


6.2845.100

Reactor plate complete to Professional Reactor

Intelligent reactor plate with reactor for PCR applications.



Phụ kiện tùy chọn

Order no.	Description	
6.1602.200	Bottle cap 2 x M6, 1 x M8 to 6.1608.120 Bottle cap for eluent bottles made of PP for eluents, post-column derivatization solutions, etc. The PP eluent bottle can withstand gas pressure when this bottle cap is used.	
6.1608.120	PP eluent bottle, 2 L Eluent bottle for applications where glass should be avoided.	
6.1836.020	Preheating capillary with core, 1.44 m Preheating capillary to Professional Reactor / Heater	
6.2061.110	Tray with sensor for Professional IC Instruments	

6.2061.120

System Connector

Connects an IC module with the footprint of a Professional IC with a Compact IC



6.2744.320

PEEK coupling M6 / UNF 10/32

Connecting transfer tubing (M6) with capillaries (UNF 10/32).



6.2845.110

HT reactor plate, complete, for Professional Reactor


Reactor plate with reactor for PCR applications at high temperatures (> 120 °C). With PEEK reactor capillaries.



6.2845.200

Reactor complete to 6.2845.100


Spare reactor for the reactor plate



6.2845.600

Column holder complete to Professional Thermostat

Intelligent column holder for the Professional Thermostat



6.2061.100

Bottle holder for Professional IC instruments

Instrument cover and bottle holder for Professional IC instruments.



6.2845.210

HT reactor, complete, for 6.2845.110

Reactor for PCR applications at high temperatures (> 120 °C). With PEEK reactor capillaries.

