



Autolab IMP

3500003360

The Autolab IMP is a compact potentiostat/galvanostat with included EIS functionality. The overall specifications of the Autolab IMP (± 10 V compliance/measured/applied voltage, ± 100 mA measured/applied current, 10 μ Hz to 1 MHz EIS frequency range and included analog integrator) make it the ideal instrument for any educational lab, as well as for starting any electrochemistry projects.

Autolab IMP is controlled by and delivered with the powerful NOVA software, including a dedicated "EC Fundamentals" procedures library, making the electrochemistry educational process even easier. The Analog and digital inputs/outputs are available to control Autolab accessories and external devices.

Phạm vi giao hàng 3500003360

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

CABLE.PWR

Power cable


Standard power cable for Autolab instruments and accessories.



1 PCS

CELLCBL.M101.204 Cell cable


Standard cell cable, 1.5 m, with connection for counter electrode (CE), reference electrode (RE), sense electrode (S), working electrode (WE) and ground for M101/M204/PGSTAT204.



1 PCS

AUT.DUMCELL.S Autolab dummy cell

Dummy cell for instrument testing.



1 PCS


CBL.USB

Standard USB cable

Standard USB cable for Autolab instruments.



Phụ kiện tùy chọn

Order no.	Description
CBL.MONIT.MAC. 204.S	Monitor cable for M101/M204/PGSTAT204 Monitor cable for M101/M204/PGSTAT204, providing connections for analog outputs (E_{out}, i_{out} and V_{out}) and input (V_{in}). 
ALL.CLIP.BLACK	Black Alligator Clamp Black alligator clamp for connections to electrodes in the electrochemical cell. 
ALL.CLIP.RED	Red Alligator Clamp Red alligator clamp for connections to electrodes in the electrochemical cell. 
AUT204.CASE	Field case for PGSTAT204 The PGSTAT204 field case is designed to conveniently store or carry the Autolab PGSTAT204 and its components. The field case for PGSTAT204 provides room for the PGSTAT204, power cable, dummy cell, cell cable, monitor cable, NOVA Getting Started Manual and NOVA installation CD. 

IME663.101.S

Interface for Metrohm 663 VA Stand

Interface for Metrohm 663 VA Stand

