



Generator electrode with diaphragm

6.00348.100

Generator electrode with diaphragm for coulometric Karl Fischer titrations, SGJ 29/22. This electrode has a cable connection angled upwards, which simplifies cable guiding and operation when used with the OMNIS Sample Robot Oven.

This generator electrode is well-suited for:

- Samples with very low absolute water content 50 ppm
- Coulometric water content determination in oil samples
- Applications where ketone reagents are used
- Applications where low-conducting reagents are used (addition of chloroform, xylene >10%)

Phụ kiện tùy chọn

Order no.	Description
-----------	-------------

6.2104.120

Electrode cable / 1 m / H

With plug H. For connecting generator electrode - KF Coulometer.



6.2104.620

Electrode cable / 0.65 m / H

With plug H, for the generator electrode - coulometer connection

