



Double Au ring electrode

6.00353.100

Glass shaft electrode with two gold rings which are polarized for redox titrations (bivoltammetric titration).

This electrode is well suited for vitamin C determination and can be easily cleaned even from very sticky samples.

Phạm vi giao hàng 6.00353.100

Qt.	Order no.	Description
10 Stk	6.1236.050	Sleeve with SGJ 14/12 mm
		Sleeve with SGJ 14/12 mm, polyethylene.



Together with 6.2008.050 storage vessel holder, it provides a support for the electrode on 807 Dosing Units.



Phụ kiện tùy chọn

Order no.	Description	
6.2104.020	Electrode cable / 1 m / F For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).	
6.02104.040	Electrode cable plug-in head G (pol.) / plug P, 0.55 m Cable for connecting an analog polarizable electrode (plug-in head G) to an analog measuring module from OMNIS.	
6.02104.050	Electrode cable plug-in head G (pol.) / plug P, 1.5 m Cable for connecting an analog polarizable electrode (plug-in head G) to an analog measuring module from OMNIS.	