



Micro Au Titrode (installation length 30.8 cm)

6.00435.120

Combined gold ring electrode with a pH glass membrane as reference electrode, installation length = 30.8 cm, narrowed to a diameter of 6.4 mm.

This maintenance-free electrode is well suited for redox titrations when the pH value remains constant, e.g.:

- determination of the chemical oxygen demand (COD, ferrometry)
- determination of penicillin and ampicillin
- for titrations using $\text{Hg}(\text{NO}_3)_2$
- for titrations in presence of chromium or iron

This electrode is stored in deion. water.

Phụ kiện tùy chọn

Order no.	Description
-----------	-------------

6.2104.020

Electrode cable / 1 m / F

For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).

