




Drain channel for OMNIS Sample Robot M/L WSM

6.01804.460

The drain channel is mounted behind the workstations of the OMNIS Sample Robot WSM M/L and, in case of an error, collects any overflowing solution in order to channel it safely and contact-free through tubing into a waste canister.

Phạm vi giao hàng 6.01804.460

Qt.	Order no.	Description
1 PCS	6.01803.000	Silicone tubing 8/12 2.5 m

Silicone tubing for connecting to the drainage of the Pick&Place module or pump module in order to channel any solution overflows that may occur safely into the waste canister.



1 PCS

6.01808.010

Y tubing connector, PP, 3 x 7-8 mm

Connector for waste tubing



1 PCS

6.1602.115

Bottle attachment / S 40 / green

For direct dosing from reagent bottles with S 40 thread (Merck...).



1 PCS

6.1624.000

Adapter SGJ 14 for 6.1619.XXX Adsorber tube

Adapter SGJ 14. Part of 6.1619.010 for Dilution Sample Processor.



1 PCS

6.2023.020

Clip for SGJ 14/15

Clip for SGJ 14/16

