



Spearhead electrode

6.0226.100

Maintenance-free pH electrode (gel electrolyte) for measurements in semisolid samples such as cheese, meat, or dough.

This electrode is stored in $c(\text{KCl}) = \text{sat.}$ (6.2308.000) and not suitable for ion-deficient solutions.

Thanks to the lifetime indicator, you always know early on when the electrode needs to be replaced.

Phạm vi giao hàng 6.0226.100

Qt.	Order no.	Description
1 PCS	6.1236.050	Sleeve with SGJ 14/12 mm
		Sleeve with SGJ 14/12 mm, polyethylene.



Together with 6.2008.050 storage vessel holder, it provides a support for the electrode on 807 Dosing Units.



Phụ kiện tùy chọn

Order no.	Description	
6.2104.020	Electrode cable / 1 m / F For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).	
6.2308.000	Electrolyte KCl (sat.), 250 mL Electrolyte KCl (sat.), 250 mL	