



## Solvotrode easyClean (fixed cable 2.0 m)

6.0229.020

Combined pH electrode with flexible easyClean diaphragm and fixed cable (2.0 m) for all automated, non-aqueous acid/base titrations.

The glass membrane is optimized for poorly conducting solutions, and thanks to the easy-to-clean easyClean diaphragm, the electrode is well suited for heavily contaminated samples, such as used oil.

This electrode can be used with non-aqueous reference electrolytes (lithium chloride or tetraethylammonium bromide). Storage in the respective reference electrolyte.

### Phạm vi giao hàng 6.0229.020

| Qt.     | Order no.  | Description                             |
|---------|------------|-----------------------------------------|
| 100 Stk | 6.1236.050 | Sleeve with SGJ 14/12 mm                |
|         |            | Sleeve with SGJ 14/12 mm, polyethylene. |



Together with 6.2008.050 storage vessel holder, it provides a support for the electrode on 807 Dosing Units.



## Phụ kiện tùy chọn

| Order no.   | Description                                                                                                                                                                                                                                                                                                                             |
|-------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 6.2320.000  | <p>TEABr 0.4 mol/L in ethylene glycol (250 mL)</p> <p>Electrolyte solution TEABr (tetraethylammonium bromide in ethylene glycol), <math>c(\text{TEABr}) = 0.4 \text{ mol/L}</math></p>                                                               |
| 6.02109.000 | <p>OMNIS adapter cable</p> <p>This adapter cable comprising socket F to plug P is suitable for use as a fixed cable (e.g. Optrode or Solvotrode easyClean) connecting a sensor without temperature sensor to the OMNIS analog measuring module.</p>  |