



Cylinder unit OMNIS special, 20 mL

6.03004.220

Intelligent 20 mL cylinder unit for one OMNIS Titrator, Titration Module or Dosing Module. This cylinder unit is especially recommended for the following solutions:

- Aqueous alkaline solutions
- Titrant 5
- Silver nitrate solutions
- Nonaqueous alkaline solutions
- Permanganate solutions
- EDTA solutions

Includes dosing tubing and antidiffusion tip.

Phạm vi giao hàng 6.03004.220

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.01508.220

Cylinder unit OMNIS special 20 mL, without accessories

Intelligent cylinder unit 20 mL for one OMNIS Titrator, Titration Module or Dosing Module. This cylinder unit is especially recommended for the following solutions:

- Aqueous alkaline solutions
- Titrant 5
- Silver nitrate solutions
- Nonaqueous alkaline solutions
- Permanganate solutions
- EDTA solutions



Without dosing tubing and antidiffusion tip.

1 PCS

6.1446.030

Ball stopper / B-14/(15)

For the movable holding of buret tips in titration vessels



2 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.1543.200

Titration tip / M6 thread

Titration tip with antidiffusion valve.



2 PCS 6.1805.100 FEP tubing / M6 / 40 cm
With light and kink protection



1 PCS 6.2008.030 Buret tip storage vessel for 807 Dosing Unit

Together with the 6.2008.050 storage vessel holder, it provides a holder for the buret tip on the 807 Dosing Unit.



1 PCS 6.2709.070 Guiding sleeve

For fixation of buret and dosing tips in titration heads and electrode holders




1 PCS 6.2739.000 Wrench

For tightening connectors



Phụ kiện tùy chọn

Order no.	Description	
6.05700.140	Consumable Kit, cylinder unit OMNIS FEP tubing connections M6 (40 cm) and antidiffusion tip for 6.03001.xxx and 6.03004.xxx OMNIS cylinder units.	
6.05700.451	Consumable Kit, cylinder unit OMNIS special 20 mL Cylinder unit special 20 mL, complete with cylinder, piston, cylinder base, valve disk, centering tube, housing, and data chip. For 6.03004.220 OMNIS cylinder unit.	