



Micro Ag Titrode

6.0433.110

Combined silver ring electrode with a pH glass membrane as reference electrode, narrowed to a diameter of 6.4 mm.

This maintenance-free electrode is well suited for precipitation titrations when the pH value remains constant (titrant silver nitrate) of, e.g.:

- chloride, bromide, iodide
- sulfides
- hydrogen sulfide
- mercaptans
- cyanides

This electrode is stored in deion. water.

Depending on the application, we recommend using an Ag Titrode with an Ag_2S coating, which can be ordered accordingly.

Phụ kiện tùy chọn

Order no.	Description	
6.2104.020	Electrode cable / 1 m / F For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).	
