



Sodium-selective glass electrode

6.0501.100

Glass electrode with sodium-selective membrane glass and SGJ.

This ISE has to be used in combination with a reference electrode and is suitable for:

- ion measurements of Na^+ (10^{-5} to 1 mol/L)
- titrations

This electrode is also suitable for difficult matrices, such as fruit/vegetable juices, and samples that contain chlorate or perchlorate.

Phụ kiện tùy chọn

Order no.	Description
-----------	-------------

6.2104.020

Electrode cable / 1 m / F

For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).



6.0750.100

LL ISE reference electrode

Silver / silver chloride reference electrode with double junction system.
This reference electrode is well suited for:

- automated applications
- ion measurements
- surfactant titrations

The ground-joint diaphragm, which is insensitive to contamination, offers a constant and reproducible electrolyte outflow. Additionally, the reference electrolyte is gelified for even better signal stability. The sensor is delivered with $c(\text{KCl}) = 3 \text{ mol/L}$ as bridge electrolyte, which can be freely selected and exchanged as needed.

