



NH₃-selective gas membrane electrode (high conc.)

Combined ammonia-selective gas membrane electrode for concentrations $>10^{-4}$ mol/L, e.g.:

- in waste water
- in samples containing oil
- for long-term measurements

An individual contaminated membrane can easily be cost-effectively replaced.

Phạm vi giao hàng

Qt.	Order no.	Description
-----	-----------	-------------

100 Stk

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



100 Stk

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Phụ kiện tùy chọn

Order no.	Description	
6.2316.030	Inner filling solution to 6.0506.100 electrode Measuring electrolyte for NH ₃ -selective electrode (50 mL)	
6.2104.020	Electrode cable / 1 m / F For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).	