



OMNIS IC-Ausr.: UF-MiPT, PP/DD

Phạm vi giao hàng

Order Number	Description	Quantity
6.02744.010	flanschlose Verschraubung M6, 1/8", 2 Stk.	2 PCS




6.01815.000	Spiralband, 3-15 mm, 0.5 m	1 PCS
-------------	----------------------------	-------



6.05000.010	Accessory kit: Bottle OMNIS IC, PP, 2 L	1 PCS
-------------	--	-------

OMNIS IC accessory set to install the reagent bottle onto the Reagent Organizer / Reagent Organizer Liquid Control.



6.2744.070

Pressure screw short

2 PCS

Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries



6.2744.040

2 x UNF 10/32 coupling

1 PCS

For connecting 1/16 in. capillaries. For IC instruments



6.2621.030

Hexagon key 4 mm

1 PCS



6.1831.060

PEEK capillary 0.5 mm i.d., 1 m

1 PCS

For inline dialysis and for Sample Changers. For IC Dialysis Unit, IC Sample Processor, Compact IC Autosampler, IC Filtration Sample Processor, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit



6.2744.034

Coupling olive/UNF 10/32, 2x

1 PCS

Connection between pressure screw and pump tubing. 2 pieces. For IC instruments with peristaltic pump.



6.1825.290

PEEK sample loop 250 µL

1 PCS

Sample loop for «Metrohm intelligent Partial loop Technique» (MiPT) and for PILS (particle into liquid sampler). With 2 Peek pressure screws.



6.2744.200

Adapter UNF 10/32 outer / M6 inner

1 PCS


Connecting 1/16 in. PEEK capillaries and transfer tubing (M6)



6.01562.100

Transferkapillare 2 mL, PEEK


1 PCS



6.01805.400

FEP-Schlauch M6, 1/8" x 1.59 mm, 3 m


1 PCS



6.01826.020

Pumpschlauch, ID 1.0 mm, weiss - weiss


1 PCS



6.01831.220

PEEK-Kap. FILTRATE OUT - INJ, 800/0.5 mm


1 PCS



6.01831.210

PEEK Kap. PP - SAMPLE IN, 450/0.5 mm

1 PCS



6.02729.110

Ultrafiltration cell 2

1 PCS



For the Inline ultrafiltration of samples together with the 858 Professional Sample Processor / 919 IC Autosampler plus. As an alternative to the time-consuming manual filtration with a syringe filter.

6.2714.020

Filtration membrane reg. cellulose, 0.2 µm, 50 pieces 1 PCS



Nominal pore size: 0.2 µm. Set of 50 pieces. For inline filtration.

6.1831.180

PEEK Capillary 0.5 mm i.d., 3 m

2 PCS



6.2744.000

Pressure screw PVDF

1 PCS



With UNF 10/32 connection. For IC instruments, inline dialysis (connection of PTFE capillaries)

6.2831.010

Plastic tweezers

1 PCS



For dialysis and filtration membranes