



OMNIS IC-Ausr.: AnCat Switch

Phạm vi giao hàng

Order Number	Description	Quantity
6.02744.010	flanschlose Verschraubung M6, 1/8", 2 Stk.	2 PCS
		
6.01806.000	PEEK capillary, 0.18 mm ID / 2 m For use in the high-pressure flow path, e.g. as preheating capillary for Metrosep Carb 2 columns in the 2 mm version (6.01090.2x0).	1 PCS
		
6.05000.010	Accessory kit: Bottle OMNIS IC, PP, 2 L OMNIS IC accessory set to install the reagent bottle onto the Reagent Organizer / Reagent Organizer Liquid Control.	1 PCS
		

6.2744.070

Pressure screw short

1 PCS


Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries



6.01831.170

PEEK-Kapillare INJ - COL, 540/0.18 mm

1 PCS



6.2744.040

2 x UNF 10/32 coupling

1 PCS

For connecting 1/16 in. capillaries. For IC instruments



6.2744.014

Pressure screw 2x

1 PCS


With UNF 10/32 connection. For the connection of PEEK capillaries



6.01805.400

FEP-Schlauch M6, 1/8" x 1.59 mm, 3 m


2 PCS



6.02845.000

Vorwärmkartusche, ID 0.18 mm

1 PCS



6.1831.180

PEEK Capillary 0.5 mm i.d., 3 m

1 PCS

