



KFT equipment for OMNIS Oven Module

6.05406.050

Complete accessory set for volumetric Karl Fischer titrations on an OMNIS Oven Module

Phạm vi giao hàng 6.05406.050

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS	6.02104.050	Electrode cable plug-in head G (pol.) / plug P, 1.5 m
-------	-------------	---

Cable for connecting an analog polarizable electrode (plug-in head G) to an analog measuring module from OMNIS.



1 PCS

6.0341.100

Double Pt-wire electrode for coulometry

Indicator electrode used for coulometric Karl Fischer titration.



1 PCS

6.1446.240

PEEK SGJ 29 stopper for KFC titration cell

Stopper for volumetric titration in a titration cell (coulometry)



2 PCS

6.1448.020

Septum 16 mm, 5 pieces

Set of 5 pcs.



1 PCS

6.1464.320

KF titration vessel / 80-250 mL / coulometric



1 PCS

6.1903.030

Stirring bar / 25 mm

Stirring bar with magnetic core, PTFE covering, length 25 mm.



1 PCS

6.2701.040

GL 18 screw cap with hole

Made of PBT. With a hole for the 6.1448.020 septum. For the 6.5405.000 cell without diaphragm and the 6.1464.32X and 6.1465.320 titration vessels.



1 PCS

6.2713.000

PTFE sleeve SGJ14

For fat-free ground-joint connection with the standard ground-joint SGJ 14.



Phụ kiện tùy chọn

Order no.	Description	
6.2811.000	Molecular sieve Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.	
6.2811.010	Molecular sieve (US) Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments. This product is intended for all users who are US residents. The material is rated hazardous in the US in accordance with OSHA Hazard Communication Standard (29 CFR 1910.1200).	