



Consumable Kit OMNIS Pump Module (2-channel)

6.05700.010

Two pump tubings each and associated press clamps for replacement in the OMNIS pump module. One Consumable Kit can be used for re-equipping two peristaltic pumps of the pump module.

Phạm vi giao hàng 6.05700.010

Qt.	Order no.	Description
2 PCS	6.01801.020	Pharm-A-line pump tubing to the peristaltic pump module

The pump tubing comes pre-equipped with adapters that permit direct installation in the peristaltic pump of the pump module. The Pharm-A-line pump tubing is used for the transport of both aqueous and nonaqueous solutions and must be operated in connection with the white press clamp in the pump module.



2 PCS

6.02703.000

Press clamp, white, for "Peristaltic" pump module

The press clamp is used to fix and secure the pump tubing in the "Peristaltic" pump module.

