



Ag/AgCl reference electrode with KCl (length 12.5 cm)

6.0726.107

Silver / silver chloride reference electrode with double junction system, installation length 10 cm. The SGJ 14/15 enables easy assembly and the flexible ground-joint diaphragm, which is insensitive to contamination, can be easily replaced. The reference and bridge electrolyte can be freely selected according to the application and are easy to exchange. This sensor is pre-filled with $c(\text{KCl}) = 3 \text{ mol/L}$ as reference and bridge electrolyte.

Phụ kiện tùy chọn

Order no.	Description
6.02105.000	Cable strand, length 0.55 m
	Electrode cables (1x banana plug 4 mm, 1x banana plug 2 mm) to connect a disconnected reference electrode to an analog measuring input in an OMNIS titration system.



6.02105.010 Cable strand, length 1.5 m

Electrode cables (1x banana plug 4 mm, 1x banana plug 2 mm) to connect a disconnected reference electrode to an analog measuring input in an OMNIS titration system, length 1.5 m.



6.1243.030 Diaphragm for reference electrode

Spare ground-joint diaphragm for 6.0726.1XX and 6.0729.1XX reference electrodes



6.2106.020 Strand / 1 m / 2 x B

Electrode cable with 2 x plug B (4 mm) for connecting electrodes with plug-in head B



6.2106.060 Strand / 2 m / 2 x B

Electrode cable with 2 x plug B (4 mm) for connecting electrodes with plug-in head B



6.2308.020

Electrolyte 3 mol/L KCl (250 mL)

Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)

