



Metrosep Organic Acids - 100/7.8

6.1005.210

The exclusion separation column for the determination of organic acids and weak mineral acids. The low capacity in comparison with the Metrosep Organic Acids - 250/7.8 (6.1005.200) allows the rapid separation of organic acids. This column is suitable primarily for small and medium concentrations in uncomplicated sample matrices.

Phạm vi giao hàng 6.1005.210

Qt.	Order no.	Description
2 PCS	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Phụ kiện tùy chọn

Order no.	Description
6.1005.250	Metrosep Organic Acids Guard/4.6 The Metrosep Organic Acids Guard/4.6 effectively removes contamination, thus protecting the analytical separation column.


