



# Metrosep C 6 - 100/4.0

6.1051.410

The 100 mm version of the Metrosep C 6 column is designed to determine standard cations, e.g. in drinking water. Excellent separation of sodium and ammonium is still achieved, despite the very short retention times. The higher capacity of the C 6 material permits larger sample volumes.

## Phạm vi giao hàng 6.1051.410

| Qt.   | Order no.  | Description      |
|-------|------------|------------------|
| 2 PCS | 6.2744.060 | Threaded stopper |

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



## Phụ kiện tùy chọn

| Order no.  | Description                                                                                                                                                                                                                                                                                                                                                                              |                                                                                     |
|------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| 6.1051.500 | Metrosep C 6 Guard/4.0                                                                                                                                                                                                                                                                                                                                                                   |  |
|            | <p>The Metrosep C 6 Guard/4.0 contains the C 6 column material and is used to protect against particles and contamination. This considerably lengthens the service life of the analytical separation column. The Metrosep C 6 Guard/4.0 works according to the "On Column Guard System" and is screwed directly onto the respective separation column with virtually no dead volume.</p> |                                                                                     |
| 6.1051.510 | Metrosep C 6 S-Guard/4.0                                                                                                                                                                                                                                                                                                                                                                 |  |
|            | <p>The Metrosep C 6 S-Guard/4.0 contains the C 6 column material and is used to protect it against particles and contamination. This considerably lengthens the service life of the analytical separation column.</p>                                                                                                                                                                    |                                                                                     |