



Metrosep Carb 2 - 150/4.0

6.1090.420

The Metrosep Carb 2 - 150/4.0 IC column is particularly suitable for the determination of carbohydrates using alkaline eluents and pulsed amperometric detection. The anion exchanger column is based on a styrene-divinylbenzene copolymer. It is stable in the range of pH = 0 to pH = 14 and separates monosaccharides and disaccharides. It is also suitable for the analysis of sugar alcohols, anhydrous sugars, oligosaccharides, etc. The column capacity has been optimized to enable the combination of rapid separations and excellent separation properties.

Phạm vi giao hàng 6.1090.420

Qt.	Order no.	Description
2 PCS	6.2744.060	Threaded stopper
For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.		



Phụ kiện tùy chọn

Order no.	Description	
6.1015.200	Metrosep BO₃ Trap 1 - 100/4.0 Trap column for the removal of borate contaminants from the eluent. The Metrosep BO₃ Trap 1 - 100/4.0 is used mainly in carbohydrate analysis with hydroxide eluents.	 A cylindrical chromatography trap column with a tan body and green end fittings. A label is visible on the side.
6.1015.300	Metrosep CO₃ Trap 1 - 100/4.0 Trap column for the removal of carbonate traces in hydroxide eluents.	 A cylindrical chromatography trap column with a tan body and green end fittings. A label is visible on the side.
6.1090.500	Metrosep Carb 2 Guard/4.0 The Metrosep Carb 2 Guard/4.0 effectively removes contamination, thus protecting the analytical separation column. The design of the guard column has been selected in such a way that its influence on the chromatographic separation can be ignored.	 A chromatography guard column with a tan body and green end fittings. It has a label that reads 'Metrosep Carb 2 Guard 4.0 100/4.0'. The design is more complex than the trap columns.
6.1090.510	Metrosep Carb 2 S-Guard/4.0 The Metrosep Carb 2 S-Guard/4.0 effectively removes contaminations, thus protecting the analytical separation column. The design of the guard column has been selected in such a way that its influence on the chromatographic separation can be ignored.	 A chromatography guard column with a tan body and green end fittings. It has a label that reads 'Metrosep Carb 2 S-Guard 4.0 100/4.0'. The design is similar to the Carb 2 Guard but with a different internal structure.