



Metrosep Carb 2 - 250/4.0

6.1090.430

The Metrosep Carb 2 - 250/4.0 IC column is particularly suitable for the determination of carbohydrates using alkaline eluents and pulsed amperometric detection. The high-capacity anion exchanger column is based on a styrene-divinylbenzene copolymer. It is stable in the range of pH = 0 - 14 and separates monosaccharides and disaccharides. It is also suitable for the analysis of sugar alcohols, anhydrous sugars, amino sugars, etc. The 250 mm version of the Metrosep Carb 2 separation column is optimized for complex separations.

Phạm vi giao hàng 6.1090.430

Qt.	Order no.	Description
2 PCS	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Phụ kiện tùy chọn

Order no.	Description	
6.1015.200	Metrosep BO3 Trap 1 - 100/4.0 Trap column for the removal of borate contaminants from the eluent. The Metrosep BO₃ Trap 1 - 100/4.0 is used mainly in carbohydrate analysis with hydroxide eluents.	
6.1015.300	Metrosep CO3 Trap 1 - 100/4.0 Trap column for the removal of carbonate traces in hydroxide eluents.	
6.1090.500	Metrosep Carb 2 Guard/4.0 The Metrosep Carb 2 Guard/4.0 effectively removes contamination, thus protecting the analytical separation column. The design of the guard column has been selected in such a way that its influence on the chromatographic separation can be ignored.	
6.1090.510	Metrosep Carb 2 S-Guard/4.0 The Metrosep Carb 2 S-Guard/4.0 effectively removes contaminations, thus protecting the analytical separation column. The design of the guard column has been selected in such a way that its influence on the chromatographic separation can be ignored.	