



Sample beaker, plastic (PP), 120 mL, 1 piece

6.1459.307

Sample beaker, plastic (polypropylene, PP), 120 mL, suitable for sample racks: 6.2041.470, 6.2041.840, 6.2041.850, 6.2068.010, 6.2068.020, 6.2068.030, and 6.9920.141, 1 piece

Phạm vi giao hàng 6.1459.307

Qt.	Order no.	Description
0 PCS	6.1459.300	Sample beaker, plastic (PP), 120 mL, 250 pieces
		Sample beaker, plastic (polypropylene, PP), 120 mL, suitable for sample racks: 6.2041.470, 6.2041.840, 6.2041.850, 6.2068.010, 6.2068.020, 6.2068.030, and 6.9920.141, 250 pieces



Phụ kiện tùy chọn

Order no.	Description	
6.1459.320	Screw cap for 61459300 / 250 pieces For safe closure of 120 mL sample beakers	
6.1459.300	Sample beaker, plastic (PP), 120 mL, 250 pieces Sample beaker, plastic (polypropylene, PP), 120 mL, suitable for sample racks: 6.2041.470, 6.2041.840, 6.2041.850, 6.2068.010, 6.2068.020, 6.2068.030, and 6.9920.141, 250 pieces	
6.2041.470	Sample rack, 22 x 120 mL, (PP) Sample rack for 810 Sample Changer, 814 USB Sample Processor, 855 Robotic Titrosampler "Light", suitable for 22 sample beakers. The following sample beakers can be used: 6.1459.300, 6.1459.307. Plastic: Polypropylene (PP)	
6.2041.840	Sample rack, 59 x 120 mL, (PVC) Sample rack for 815 Robotic USB Sample Processor XL, 855 Robotic Titrosampler "Basic", suitable for 59 sample beakers. The following sample beakers can be used: 6.1459.300, 6.1459.307. Plastic: Polyvinyl chloride (PVC)	

6.2041.850 Sample rack, 59 x 120 mL, (PP)

Sample rack for 815 Robotic USB Sample Processor XL, 855 Robotic Titrosampler "Basic", suitable for 59 sample beakers. The following sample beakers can be used: 6.1459.300, 6.1459.307.


Plastic: Polyvinyl chloride (PVC)



6.2068.010 Sample rack, 20 x 120 mL, (PP)

Sample rack for 864 Balance Sample Processor, suitable for 20 sample beakers. The following sample beakers can be used: 6.1459.300.

Plastic: Polypropylene (PP)



6.2068.020 Sample rack, 24 x 120 mL, (PVC)

Sample rack for 815 Robotic Soliprep, suitable for 24 sample beakers. The following sample beakers can be used: 6.1459.300, 6.1459.307.

Plastic: Polyvinyl chloride (PVC)



6.2068.030 Sample rack, 10 x 250 mL, (PVC)

Sample rack for 815 Robotic Soliprep, suitable for 10 sample beakers. The following sample beakers can be used: 6.1432.320, 6.1432.323, 6.1453.220, 6.1453.250.



6.9920.141

Sample rack, 53 x 120 mL, (PP)

Sample rack for 815 Robotic USB Sample Processor XL, 855 Robotic Titrosampler "Basic", suitable for 53 sample beakers. The following sample beakers can be used: 6.1459.300, 6.1459.307.

Plastic: Polypropylene (PP)



6.1458.040

Titration Head, 3 x SGJ 14

Titration head for 120 mL beakers. With 3 openings SGJ 14 and 4 openings 6.4 mm, 2 dosing tips included. Suitable for rinsing and aspiration equipment.



6.2037.080

Foil holder

Foil holder for 120 mL PP beaker (6.1459.300)



6.1909.050

Stirring propeller / 102 mm

Propeller for 802 Stirrer when this is used with sample racks for 120 mL beakers (6.1459.300)



6.1909.060

Stirring propeller / 95 mm

Stirring propeller with 3 wings for intensive stirring

